

Bridging Signal and Intelligence: A DSP-Driven EEG Feature Pipeline for Neural Network Seizure Detection

Presented By

Md. Ohiduzzaman Sajol (2003043)
Durjoy Sarkar Dhrubo (2003026)

Presented To

Dr. Monira Islam
Associate Professor, Dept. of EEE, KUET

Presented On

November 2, 2025



Outline

- **Abstract**
- **Introduction**
- **Background**
- **Methodology**
- **Results**
- **Discussion**
- **Conclusion**
- **Future Work**
- **References**

Abstract

- **DSP-based EEG pipeline for seizure detection**
- **CHB-MIT dataset with notch & bandpass filters**
- **Extracted spectral & statistical EEG features**
- **Visualized raw vs. filtered & frequency spectra**
- **Features ready for ML seizure classification**

Introduction

- EEG is a **non-invasive** method to measure brain activity.
- Automated EEG-based **seizure detection** is a key biomedical DSP task.
- EEG's **non-stationary** nature demands **strong preprocessing & feature extraction**.
- DSP provides tools for **time–frequency** domain analysis (e.g., **convolution, DFT, filtering**).
- EEG data (**CHB-MIT**) is preprocessed, features extracted, and prepared for ML seizure classification.

Background: EEG and Seizure Activity

- EEG records brain's electrical signals via scalp electrodes.
- Typical EEG range: **0.5–70 Hz**, divided into standard bands:
 - Δ (**0.5–4 Hz**) – Deep sleep
 - θ (**4–8 Hz**) – Drowsiness/meditation
 - α (**8–13 Hz**) – Relaxed wakefulness
 - β (**13–30 Hz**) – Active focus
 - γ (**30–70 Hz**) – Cognitive processing
- **Seizures:** sudden, synchronized bursts of activity
 - visible as amplitude spikes & spectral power changes

Methodology: Data Acquisition

- Dataset: **CHB-MIT Scalp EEG Database** (pediatric epilepsy patients)
- Format: **EDF files** with multi-channel EEG recordings
- Loaded via *mne.io.read_raw_edf()*
- Includes both **seizure** and **non-seizure** intervals

Methodology: Preprocessing Pipeline

- **Notch Filter (60 Hz):** remove power-line noise (IIR notch)
- **Bandpass FIR Filter (0.5–70 Hz):** isolate standard EEG range
- **Segmentation:** 10-s overlapping epochs for temporal analysis
- **Normalization:** z-score (zero mean, unit variance) for consistency

Methodology: DSP Operations

- **Correlation:** inter-channel Pearson correlation matrices
- **Convolution:** moving-average smoothing to reduce fluctuations
- **Spectral Analysis:** FFT + Welch PSD for frequency and band-power estimation

Methodology: Feature Extraction

- **Band Powers:** mean, std, median in δ , θ , α , β , γ bands
- **Statistical:** energy, variance, skewness, kurtosis
- **Spectral Entropy:** Shannon entropy of normalized PSD
- **Correlation Metrics:** mean, std, max of inter-channel correlations
- **Band Ratios:** e.g., α / β for cognitive state insight

Methodology: Implementation

- Tools: **Python 3.10**, MNE, SciPy, NumPy, pandas, Matplotlib
- Interactive workflow via **Jupyter Notebook**
- Features aggregated into *chb_features.csv*
- Repository: github.com/projectohid/4142-project

Results: Spectrogram Analysis (Normal vs. Seizure)

- **Normal EEG (Fig. 1)** : Broad spectral power (10–30 Hz, α & β bands)
- **Seizure EEG (Fig. 2)** : Dominant low-frequency power (<10 Hz, δ & θ)
- Confirms effective filtering and **distinct physiological vs. pathological patterns**

Results: Spectrogram Analysis (Normal vs. Seizure)




Fig. 1: Spectrogram for *Chb01_01.edf* (non-seizure).

Results: Spectrogram Analysis (Normal vs. Seizure)




Fig. 2: Spectrogram for *Chb01_03.edf* (seizure).

Results: Raw Vs. Filtered Signal

- **Time-domain (Fig. 3):** Filtered signal smoother; reduced drift & 60 Hz noise
- **Frequency-domain (Fig. 3):** 60 Hz peak suppressed; spectrum confined to 0.5–70 Hz
- Verifies **notch + FIR bandpass filters** effectiveness

Results: Raw Vs. Filtered Signal




Fig. 3: Raw and filtered Signals and their Spectrums for *Chb01_01.edf*.

Results: EEG Band Power Distribution

- Power hierarchy:

$$\Delta > \theta > \beta > \gamma > \alpha$$

- Low-frequency dominance
→ seizure or drowsy state indicator
- Confirms **spectral accuracy**
of feature extraction




Fig. 4: Average Band Power Across EEG Frequency Bands for *Chb01_01.edf* .

Results: Time-Domain Statistical Features

- **High kurtosis:** sharp, transient seizure spikes
- **Low skewness:** symmetric yet with frequent outliers
- Confirms **impulsive, non-Gaussian EEG** patterns in seizures




Fig. 5: Distribution of Statistical Time-Domain Features for *Chb01_01.edf*.

Results: Spectral Entropy vs. Delta Power

- **Negative correlation** → higher δ power = lower entropy
- Indicates **reduced complexity** during seizure activity
- Validates **entropy as nonlinear brain state marker**




Fig. 6: Spectral Entropy vs. Delta Band Power for *Chb01_01.edf*.

Results: Inter-Channel Correlation

- **Max correlation ≈ 1.0 :**
strong cross-channel
synchronization
- **Consistent** mean correlation
= stable phase coupling
- Demonstrates **spatial
coherence** in seizure EEG




Fig. 7: Distribution of Inter-Channel Correlation Metrics for *Chb01_01.edf*.

Results: Feature Matrix & ML Readiness

Table I: Part of the extracted Feature Matrix for *Chb01_01.edf*.

	delta_mean	delta_std	delta_max	delta_median	theta_mean	theta_std	theta_max	theta_median	alpha_mean	alpha_std
0	0.439917	0.122878	0.680715	0.425868	0.134942	0.048450	0.229849	0.132495	0.046567	0.029513
1	0.491996	0.147472	0.753315	0.489125	0.157016	0.071678	0.355560	0.139707	0.052553	0.036122
2	0.375810	0.082527	0.485999	0.389008	0.155521	0.067377	0.320863	0.145871	0.055865	0.031640
3	0.328070	0.089442	0.475458	0.338504	0.164152	0.065626	0.293830	0.164945	0.056487	0.021148
4	0.467772	0.085527	0.657703	0.458547	0.173338	0.052196	0.284436	0.168564	0.054683	0.029469

5 rows × 34 columns

- Each 10 s epoch → **30+ DSP-based features**
- Seizure/non-seizure epochs clearly separable in feature space
- Confirms **pipeline suitability for neural network models**

Results: Summary

Domain	Key Observation	DSP Technique
Spectrograms	Low-frequency dominance	STFT / Spectrogram
Filtering	Removed 60 Hz & artifacts	Notch + FIR
Band Power	High δ , θ ; low α , β , γ	FFT + PSD
Statistics	High kurtosis, low skewness	Time-domain
Entropy	–ve correlation with δ power	Entropy calc.
Correlation	Strong inter-channel sync	Pearson corr.
Feature Set	Clear seizure separation	DSP feature extraction

Discussion

- **Spectral Behavior**

- Normal EEG: broad power (10–30 Hz, α – β bands)
- Seizure EEG: strong low-freq. dominance (<10 Hz, δ – θ)
- Confirms **shift from desynchronized** → **synchronized** neuronal activity
- Filtering (notch + bandpass) effectively removed 60 Hz noise & drift while preserving core EEG content

- **Band Power Trends**

- Hierarchy: $\Delta > \theta > \beta > \gamma > \alpha$
- High α/β ratio & reduced high-freq. power → decreased cortical processing during seizures

Discussion

- **Time-Domain Statistics**
 - **High kurtosis:** sharp epileptic transients
 - **Low skewness:** symmetric but heavy-tailed amplitude distribution
 - Indicates **impulsive, non-Gaussian** seizure morphology
- **Entropy & Complexity**
 - **Entropy** \downarrow as **Δ -power** $\uparrow \rightarrow$ reduced signal complexity during seizures
 - Captures transition to **predictable, synchronized** brain dynamics
 - Entropy = nonlinear complement to power-based metrics

Discussion

- **Inter-Channel Correlation**
 - **High cross-channel correlation (≈ 1.0)** → strong spatial coherence
 - Reflects widespread cortical synchronization in seizure states
- **Integrated DSP Takeways**
 - **Filtering:** ensures clean, physiologically valid signals
 - **FFT / PSD:** reveal spectral composition
 - **Statistics + Entropy:** quantify morphology & complexity
 - **Correlation:** captures spatial coupling → Together, **classical DSP techniques uncover physiological + pathological EEG patterns**, producing interpretable features for intelligent seizure detection.

Conclusion

- DSP methods effectively applied for EEG seizure analysis.
- Filtering + spectral + statistical tools made EEG **interpretable**.
- **Key Points:**
 - **Preprocessing:** Notch & FIR filters removed noise, kept key signals.
 - **Spectral:** Clear normal–seizure separation (low-freq dominance).
 - **Stats:** High kurtosis, low entropy → synchronized seizure waves.
 - **Correlation:** Strong inter-channel coupling during seizures.
 - **Features:** Compact, ML-ready dataset for classification.
- **Takeaway:**
 - DSP bridges theory & application, enabling explainable EEG systems.

Future Work

- **Wavelet Filters:** Handle dynamic, non-stationary EEG noise.
- **ICA Cleaning:** Remove eye/muscle artifacts.
- **Deep Learning:** Train models on DSP-derived features.
- **Automation:** Streamline preprocessing & feature selection.

References

- [1] L. Tan and J. Jiang, *Digital Signal Processing: Fundamentals and Applications*, 3rd ed. Elsevier Science, 2018.
- [2] M. Teplan, “Fundamentals of EEG measurement,” *Measurement Science Review*, vol. 2, no. 2, pp. 1–11, 2002.
- [3] H. Adeli, Z. Zhou, and N. Dadmehr, “Analysis of EEG records in an epileptic patient using wavelet transform,” *Journal of Neuroscience Methods*, vol. 123, no. 1, pp. 69–87, 2003.
- [4] U. R. Acharya, S. V. Sree, G. Swapna, R. J. Martis, and J. S. Suri, “Automated EEG analysis of epilepsy: A review,” *Knowledge-Based Systems*, vol. 45, pp. 147–165, 2013.
- [5] A. Aarabi, R. Fazel-Rezai, and Y. Aghakhani, “EEG seizure prediction: Measures and challenges,” *Computers in Biology and Medicine*, vol. 42, no. 4, pp. 376–390, 2012.

Thank You

Any Questions?

